


D-[1-¹³C]arabinose
¹³C {¹H} Spectrum (75 MHz) in D₂O





Assignments

α-pyranose
C1 98.2
C2 73.4
C3 74.0
C4 70.0
C5 67.9

β-pyranose
C1 94.1
C2 70.0
C3 70.1*
C4 70.2*
C5 64.0

α-furanose
C1 102.6

β-furanose
C1 96.6



Specialty Stable Isotopes for the Biomedical Sciences

115 South Hill Street
South Bend, Indiana, USA 46617-2701
www.omicronbio.com

Phone: 574-287-6910
Fax: 574-287-7165
Email: omicron@omicronbio.com