



D-[1-¹³C]lyxose
¹³C {¹H} Spectrum (75 MHz) in D₂O



Assignments

α -pyranose	
C1	95.5
C2	71.5
C3	72.0
C4	69.0
C5	64.6

β -pyranose	
C1	95.9
C2	71.5
C3	74.2
C4	68.0
C5	65.7



Specialty Stable Isotopes for the Biomedical Sciences

115 South Hill Street
South Bend, Indiana, USA 46617-2701
www.omicronbio.com

Phone: 574-287-6910
Fax: 574-287-7165
Email: omicron@omicronbio.com